


Optimizing CMN collection for ECP With OPTIA device by filling one single chamber



SFH Paris December 12th 2014
Dr Kanouni Tarik
Therapeutic Apheresis Unit
University Hospital Saint Eloi Montpellier



Optimizing CMN collection for ECP With OPTIA device by filling one single chamber

⊕ Rational of the study

⊕ Large variability of CMN collected

- Many patient had few CMN collected
- Specialty those with immunosuppressive therapy (cortison with low level of CMN in the blood
- CMN collected depending of CMN in the blood

⊕ No target for amount of CMN or lymphocyte

⊕ What is the current procedure.

- To process 1 or 2 total blood volume wirh cells separator
- Therakos 1.5 liter processed


⊕ Aim of this study

⊕ Based on the fact that we don't know if total amount of CMN have an impact of therapy,

⊕ Not to collectec as much we can

⊕ But to Standardize CMN amount for each ECP session

- All patient on the same level regarding on CMN collection
- For those patient to hom we collete few CMN cells
- Easy to make any retrospectiv analysis on clinical ECP clinical effect, having one less variability data.



Our experience on CMN collection in our unit over more than 1800 ECP session

1/ Wide experience on device used

Cobe Spectra



OPTIA



Amicus




COMTEC

Stemle connection done before CMN collection


UVA PIT System



Thérakos CellEx



UVAR XTS



2/ Exhaustive quality control on CMN

CBC and differential

IN the blood before ECP

In collection bag

Real time integrated in Apheresis database







Epic

UVA PIT Spectra


-25-



1/ Large variability of CMN collected



2/ CMN collected depending of CMN in the blood




How to reach the goal of the study

✚ Looking for technique:


- ✚ Easy , repeatable


- ✚ Fully automated (reduce human variability)

- ✚ By using OPTIA Device we can reach this target




Why using OPTIA with one chamber






Quality of cells trapped in chamber



- ⊕ WBC in chamber:
 - ⊞ CMN: Average is 91%
 - ⊞ Mediane: 98%
- ⊕ Almost all WBC in the chamber are CMN
- ⊕ When a chamber is filled
 - ⊞ Should we have always the same amount of cells?



Study parameter

1 single chamber setting

- Set UP One single chamber in the software of OPTIA

- 40 ml in total in the bag

- 16 ml as collected chamber volume
- 24 ml to rinse the tubing so that we get enough volume for processing into PIT systeme

From 01/2013 to 04/2014

- 146 ECP Sessions

- 102 complete for analysis

Comparison:

- 202 ECP done by Optia


- TBV between 0.9 and 1.1 processed

Data to be compared

- CMN, Lymp, Mono


- Variability

- Duration, Total blood volume preprocessed,




Results: Regarding Total CMN collection

- ⊕ Regarding to variability
 - ⊞ Reducing by
flilong on single
chamber
- ⊕ Regarding to total
cell collected
 - ⊞ Effective




Intra individual variation

- ✚ Analysis for intraindividual patient variability
 - ▣ Each point is patient
 - ▣ Sorted by number of procedure
 - ▣ Analysis for total CMN collection
- ✚ By filling 1 single chamber
 - ▣ Reduce intraindividual variation
 - ▣ Variability is more predictable




*Analysis of the amount of total patient
blood volume processed*



Duration of CMN collection



	8	9	10	11	12	13	14	15	16	17	18
'Salle 1'											
Chb:01			<div style="border: 1px solid black; padding: 2px;"> RA MI 08:53>09:46 </div>	<div style="border: 1px solid black; padding: 2px;"> UV 33,5 Min </div>			<div style="border: 1px solid black; padding: 2px;"> RU PH 12:14>13:25 (D=01:11) </div>	<div style="border: 1px solid black; padding: 2px;"> UV 24 Min </div>			
'Salle 2'											
Chb:02			<div style="border: 1px solid black; padding: 2px;"> BO Ab 09:27>10:21 </div>						<div style="border: 1px solid black; padding: 2px;"> BA AN 14:03>15:55 (D=01:52) ##CELLE## </div>	<div style="border: 1px solid black; padding: 2px;"> UV 25 Min </div>	
'Salle 3'											
Chb:03				<div style="border: 1px solid black; padding: 2px;"> BE SA 10:10>12:11 (D=02:01) ##COMTE## </div>		<div style="border: 1px solid black; padding: 2px;"> MA FR 13:25>15:17 (D=01:52) ##COMTE## </div>					
'Salle 4'											
Chb:04			<div style="border: 1px solid black; padding: 2px;"> CA MA 09:50>10:43 </div>				<div style="border: 1px solid black; padding: 2px;"> JO Ch 13:15>15:39 (D=02:24) </div>				
'Salle X'											
Chb:05											
<div style="display: flex; justify-content: space-between;"> LUNDI 24/03-(9) MARDI 25/03-(8) MERCREDI 26/03-(10) JEUDI 27/03-(9) VENDREDI 28/03-(9) SAMEDI 29/03-(1) DIMANCHE 30/03-(0) </div>											



Conclusion

✚ Filling 1 single chamber as target

- ▣ Is effective way according to standardize CMN collection
- ▣ Easy to manage (always the same volume collected)
- ▣ Good for schedule management shorter and predictable
- ▣ Getting enough cells

✚ To collect more CMN cells

- ▣ To maintain the advantages of chamber target.
- ▣ We can increase the target to 2 chambers
- ▣ Knowing that the second chamber is filled faster than the first

✚ This is intermediate analysis

- ▣ We have to check the impact of
 - Platelets, RBC cells, individual
- ▣ Waiting for technical analysis from Terumo
 - Each session data was sent to Terumo R&D to check if there was any technical failure during collection to make this data more reliable